
Panasonic

Portable Oral Irrigator
Irrigador oral portátil
Irrigateur bucal portable

Model No. EW176
Modelo No. EW176
No. de modèle EW176

Operating Instructions Instrucciones de funcionamiento Mode d'emploi

Thank you for purchasing the Panasonic Portable Oral Irrigator EW176.
Muchas gracias por adquirir el irrigador oral portátil Panasonic EW176.
Nous vous remercions d'avoir acheté l'irrigateur bucal portable Panasonic modèle EW176.
Before operating this device, please read these instructions completely.
Antes de usar este dispositivo, lee completamente estas instrucciones.
Avant de l'utiliser, veuillez lire attentivement tout ce mode d'emploi.






INTRODUCTION

Thank you for purchasing the Panasonic Oral Irrigator EW176. With its water jet of 1,600 pulses per minute, the Panasonic Oral Irrigator gently yet thoroughly cleanses your teeth and gums, removing trapped food particles and reducing bacteria that can lead to gingivitis and other periodontal diseases. When used in combination with regular brushing and flossing, the Panasonic Oral Irrigator can help you have healthier teeth and gums.


Before using this device, please read these instructions completely. Keep these instructions for reference.

IMPORTANT SAFETY PRECAUTIONS


1. Do not use while bathing.
2. Do not place the device where it can slide, fall or be pulled into water.
3. Close supervision is necessary when this device is used by, on or near children.
4. **KEEP JET NOZZLE ATTACHMENT AWAY FROM SMALL CHILDREN.** The nozzle is small enough to be swallowed and can become lodged in the throat.
5. Keep this device out of the reach of children when not in use.
6. This device is designed for use only in the oral cavity. Do not direct water in delicate areas, such as under the tongue, in the eye, ear or nose, or in other delicate areas of the body.
7. Use only attachments recommended by the manufacturer.
8. **DO NOT** operate this product if it is not working properly or it has been dropped or damaged. If damaged or inoperative, return the product to a service center for inspection and repair.
9. **NEVER** drop or insert any object into any opening on the device.
10. Do not disassemble the device. Damage to the product and/or personal injury may occur.
11. People with tooth or gum disease or people receiving oral treatment should consult their dentist before use.
12. People who experience tooth or gum pain may not be able to use this device and should consult their dentist.
13. **NEVER** use while drowsy.
14. Remove the batteries from the device if it will not be used for more than a month. Otherwise, electrical leakage may occur if the insulator is worn.
15. Fill the reservoir with water only.
16. Certain types of mouthwash may damage the equipment by causing the housing or water reservoir to crack. Therefore, do not use mouthwash in the water reservoir.


6

7

8

Index

ENGLISH	2
ESPAÑOL	10
FRANÇAIS	16

PARTS IDENTIFICATION

(Fig. 1)

Main unit

- ① Nozzle insertion hole
- ② Cap
- ③ Battery cover
- ④ Nozzle release button
- ⑤ Lock button
- ⑥ Switch
- ⑦ Nozzle compartment
- ⑧ Water reservoir

- There may be moisture inside the water reservoir and suction hose. This is moisture left from the distilled water used for product testing at the factory and is hygienically safe.

- ⑨ Water reservoir lid

Accessory

- ⑩ Nozzle (1)

BEFORE USE

Loading Batteries

- Be sure to wipe the unit before installing the batteries so that no water enters the battery compartment.
1. Align the cap with the mark on the unit as shown in Fig. 2 and remove the cap.
 2. Twist the battery cover to the OPEN position and open it (Fig. 3).
 3. Insert two 1.5-V alkaline (AA/LR6) batteries, ensuring correct "+" and "-" polarity.
 - The device can be used approximately 60 times (with one full reservoir of water per use) with a fresh set of Panasonic alkaline batteries.
 - With rechargeable batteries (Panasonic Super Ni-Cd 1000/DC 1.2-V, P-3SPS) installed, the device can be used approximately 60 times.

- Do not mix old and new batteries.

4. Twist the battery cover to the CLOSE position and close it.
5. Replace the cap, aligning the mark on the unit as shown in Fig. 4.

Note:

The batteries' power output is impaired at low temperatures.

CAUTION:

Do not throw batteries into a fire. Doing so may result in an accident or explosion.

USING YOUR ORAL IRRIGATOR (Figs. 5 & 6)

- Rinsing your mouth with mouthwash is not a substitute for proper dental care. Brush your teeth or dental prosthesis for at least five minutes immediately after every meal, then use the device.
1. Turn the cap to align its hole with the unit's nozzle insertion hole.
 2. Remove the nozzle from its compartment and insert firmly into the nozzle insertion hole until the nozzle clicks into place.
 3. Twist the water reservoir section to the left and pull it downward to the water reservoir usable line.
 4. Open the water reservoir lid and fill the reservoir with cold or warm water, not warmer than 104°F (40°C).
 - Do not put any toothpaste, powder or other chemicals into the water reservoir.
 - Certain types of mouthwash may damage the unit by causing the housing or reservoir to crack. Therefore, use only water in the water reservoir.

CAUTION

NEVER fill the water reservoir with hot water! Sensitive gums can be easily irritated or damaged by hot water.

5. Close the lid.
6. Lean over the sink and place the nozzle in your mouth. Close your lips slightly to prevent splattering.
 - Hold the device upright when in use. Excessive tilting can prevent the water from coming out of the nozzle.
7. Pressing the lock button, slide the switch up to turn on the device.
 - Be careful not to press the nozzle release button during use. If the release button is pressed, the water pressure may push and release the nozzle.

TIPS ON USE

1. Bend over the sink and hold the unit in an upright position with the nozzle tip pointing toward your teeth.
2. Keep your mouth opened slightly so that the water can flow out.
3. Starting in the back (molar) area and working toward the front teeth, systematically rinse all teeth and interdental spaces by:
 - a) directing the water jet flow at a right angle to the teeth and the area between the teeth and gums;
 - b) directing the water jet slowly along the teeth;
 - c) directing the water jet along the gumline with the nozzle tip at a right angle to the gum tissue;
 - d) cleaning between the teeth and around any braces, crowns or bridgework thoroughly.
 - In addition to using the oral irrigator and brushing and flossing your teeth, we recommend regular checkups by a dentist as part of a complete dental hygiene program.

PRECAUTION

- Do not turn the switch on before putting the nozzle in your mouth.
- Do not press the tip of the nozzle too firmly against your teeth or gums. Doing so may injure your gums.
- If you are suffering from tooth or gum disease or receiving oral treatment, consult a dentist before use.
- The device is designed exclusively for use in the oral cavity. Never use it for any other purpose.
- During the first several uses, light bleeding of the gums may occur. If your gums are healthy, any bleeding should stop in one to two weeks. If bleeding persists for two or more weeks, consult your dentist before continuing use.

AFTER USE (Fig. 7)

1. Open the water reservoir lid and pour out the water.

Caution: NEVER immerse the device in water.
2. Turn on the device for a few seconds until any water remaining inside the device and nozzle drains out completely.
3. Close the water reservoir lid.
4. Push the reservoir up and twist it to the right to lock it in place.
 - Be careful not to accidentally insert or pinch your finger.
5. Pressing the nozzle release button, remove the nozzle and place it into the nozzle compartment.
6. Turn the cap and secure the nozzle.
 - After rinsing the unit, wipe and dry it well before storing the unit to prevent contamination and bacterial growth.

MAINTENANCE

① Nozzle

Remove the nozzle from the device by pressing the nozzle release button. Rinse with water and wipe with a dry cloth.

② Cap

Remove the cap from the device by aligning the mark as shown in Fig. 8-a. Rinse with water and wipe with a dry cloth.

③ Main unit

Clean the main unit by wiping it with a cloth lightly dampened with liquid soap.

- Use only water or a neutral cleaning agent to clean the main unit. Never use caustic or abrasive agents (e.g., vinegar or lime-scale remover) because they may damage the device.

④ Reservoir

Align the reservoir replacing marks as shown in Fig. 8-b and pull the reservoir down to remove it.

Clean the inside with a soft brush or similar object if necessary. Wipe off any residue inside the reservoir. If the device will not be used for more than a week, be sure to dry the inside of the reservoir.

⑤ Water reservoir lid

Wipe off to prevent water leakage.

⑥ Hose

Be careful not to bend the hose. Doing so may damage the hose and result in a malfunction.

CAUTION

- Do not wash with hot water (exceeding 122°F or 50°C).
- Do not immerse the main unit in water.

Note:

If the device will not be used for more than a month, be sure to remove the batteries.

TRAVELING WITH THE PANASONIC ORAL IRRIGATOR

Thanks to its handy shape and compact size, the Panasonic Oral Irrigator is an ideal traveling companion. When using it on a trip, be sure to completely drain the reservoir and nozzle before putting the device away.

TROUBLESHOOTING -- Possible Causes and Remedies

Problem	Possible cause	Remedy
The device stops immediately after being switched on or does not operate at all.	Batteries are not inserted properly.	Insert the batteries properly. See page 6.
	Battery power is low.	Replace the batteries with new ones. See page 6.
Water pressure is too low or operating time is very short.	Battery power is low.	Replace the batteries with new ones. See page 6.
The pump is running but no water comes out of the nozzle.	The reservoir is empty.	Fill up reservoir with water.
	The reservoir is not pulled down enough.	Pull down the reservoir fully.
	Wrong position.	Use the device in an upright position with the nozzle pointing upwards.
	The nozzle is clogged or deformed.	Replace the nozzle.

If the device still does not function properly, stop using and contact an authorized service center or the store from which the device was purchased.

CAUTION

- Never attempt to repair the device yourself.
- If parts must be exchanged or replaced, or if the device malfunctions, contact an authorized service center.

SPECIFICATIONS

Power supply:	3 V DC (Two 1.5-V alkaline batteries (AA/LR6))
Pump frequency:	1600 pulses/minute
Water pressure:	approx. 440 kPa (approx. 4.5 kgf/cm ²)
Battery life:	approx. 20 minutes (with Panasonic alkaline AA/LR6 batteries)
Reservoir volume:	90 ml (approx. 20 seconds of use with new batteries)
Dimensions:	7-3/4" (H) x 1-3/4" (diameter.) (194 (H) x 46 (diameter) mm) 5-1/4" (H) x 1-3/4" (diameter.) (135 (H) x 46 (diameter) mm) when the tank is closed.
Weight:	6 oz (170 g) excluding batteries

Specifications are subject to change without notice.

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Matsushita Electric Works, Ltd.